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Print...

Print This Page

[Important Medical Record Information](#)

ICD Model: Marquis DR 7274

Clinical Status: Since Dec 10, 2000

Episodes		% Pacing	
VF	0	AS-VS	<0.1 %
FVT	0	AS-VP	0.0 %
VT	0	AP-VS	99.9 %
SVT	0	AP-VP	0.1 %
NST and Others	0		
Mode Switch	0		

Observations (0)
- No observations based on currently interrogated data.

ICD Status
Battery Voltage (ERI=2.62V) 3.1
Last Full Energy Charge 5.29
Last Capacitor Formation (Interval=Auto)

Lead Information

Atrial	Medtronic	5076 CapSureFix(R) Novus	PJN1234567
Ventricular	Medtronic	6947 Sprint Quattro Secure (TM)	TDG123456H

Lead Performance	Atrial	Ventric
EGM Amplitude	2.4 mV	15.0 mV
Pacing Impedance	480 ohms	776 ohms
Defibrillation Impedance		56 ohms
SVC (HVX) Impedance	64 ohms	

Adobe Acrobat Professional - [Concerto-All.pdf]

File Edit View Document Comments Tools Advanced Window Help

Search Create PDF Comment & Markup Send for Review Secure Sign Forms

Select 80% Help

Bookmarks

- Smith, John. Concerto™ C154DW
 - Current EGM
 - Quick Look
 - Cardiac Compass™ with Non-su
 - Cardiac Compass™ with Heart F
 - Battery and Lead Status
 - Lead Performance Trends
 - Counters
 - Parameters
 - Arrhythmia Episode List
 - Treated VT/VF Episode #2
 - episodePlot
 - episodeEGM
 - episodeText
 - Non-sustained VT Episode #6
 - episodePlot
 - episodeEGM
 - episodeText
 - V. Sensing Episode List
 - Patient Alert
 - Rate Histograms
 - Patient Information
 - Heart Failure Management

Signatures

Pages

Attachments

Comments

Medtronic

Quick Look

Device: Concerto™ C154DWK Serial Number: PVR100012Q Date of Interrogation: 26-Jun-2006 00:27:18
Patient: Smith, John Physician: Kelly, Sarah

Device Status

				Measured on:
Battery Voltage (RRT=2.62 V)	3.20 V			26-Jun-2006
Last Full Charge	7.29 sec			02-Jun-2006

Pacing Impedance

	Atrial	RV SVC	LV	
Pacing Impedance	472 ohms	480 ohms	376 ohms	02-Jun-2006
Defibrillation Impedance		RV=51 ohms SVC=63 ohms		02-Jun-2006

Capture Threshold

Programmed Amplitude/Pulse Width	3 V / 0.4 ms	3 V / 0.4 ms	1 V at 0.4 ms 4 V / 0.4 ms	

Measured P/R Wave

Measured P/R Wave	2.3 mV	22.4 mV		02-Jun-2006
Programmed Sensitivity	0.45 mV	0.3 mV		

Parameter Summary

Mode	DDDR	Lower Rate	80 bpm	Paced AV	110 ms
Mode Switch	150 bpm	Upper Track	130 bpm	Sensed AV	100 ms
V. Pacing	RV->LV	Upper Sensor	120 bpm		

Detection

		Rates	Therapies
AT/AF	Monitor	>150 bpm	
VF	On	>188 bpm	ATP During Charging. 35J x 6
FVT	via VF	188-214 bpm	Burst(1), 35J x 5
VT	On	150-188 bpm	Burst(3), 5J, 25J, 35J x 3

Enhancements On: AF/Afl, Sinus Tach, 1:1 SVT

Clinical Status


	Since 02-Jun-2006	Cardiac Compass Trends (May-2006 to May-2007)
Treated		
VF	0	
FVT	0	
VT	1	
Monitored		
VT (Off)		
VT-NS (>4 beats, >150 bpm)	5	
SVT: VT/VF Rx Withheld	2	
AT/AF	0	
Time in AT/AF	<0.1 hr/day (<0.1%)	
Functional		
Patient Activity	Last Week <2 hr/day	

Cardiac Compass Trends (May-2006 to May-2007)

Treated ≥5
VT/VF (#/day)
AT/AF (min/day)
Patient Activity (hr/day)

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
Adobe Acrobat Professional - [Concerto-All.pdf]

File Edit View Document Comments Tools Advanced Window Help

Create PDF Comment & Markup Send for Review Secure Sign Forms

Select 80% Help

Bookmarks

- Smith, John. Concerto™ C154DW
 - Current EGM
 - Quick Look
 - Cardiac Compass™ with Non-su
 - Cardiac Compass™ with Heart F
 - Battery and Lead Status
 - Lead Performance Trends
 - Counters
 - Parameters
 - Arrhythmia Episode List
 - Treated VT/VF Episode #2
 - episodePlot
 - episodeEGM
 - episodeText
 - Non-sustained VT Episode #6
 - episodePlot
 - episodeEGM
 - episodeText
 - V. Sensing Episode List
 - Patient Alert
 - Rate Histograms
 - Patient Information
 - Heart Failure Management

Signatures

Pages

Attachments

Comments

Medtronic

Treated VT/VF Episode #2

Device: Concerto™ C154DWK Serial Number: PVR100012Q Date of Interrogation: 26-Jun-2006 00:27:18
Patient: Smith, John Physician: Kelly, Sarah

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












		LATITUDE SM Patient Management				HELP/CONTACT US		LOG OUT	
		Patients For Review Hood, Maynard F				Logged In: Hood, Maynard F 20 Jan 2005			
PATIENTS FOR REVIEW		ALL PATIENTS		CONFIGURE					
 Print Reports  Review Complete  Dismiss Patient									
☐ All	Patient/Device	Last Remote Interrogation	My Alerts	Device Status	Physician(s)	Sched.	Disposition		
☐	Becker, Alton RENEWAL 3 RF H215/539610	28 Sep 2004			Hood, Maynard F	Yes	Review Started		
☐	Ecker, Ute H RENEWAL 3 RF H215/535755	28 Sep 2004			Hood, Maynard F	Yes	New Patient New PG		
☐	Gordon, Debbie RENEWAL 3 RF H215/538706	28 Sep 2004			Hood, Maynard F	Yes	New Patient New PG		
☐	Lippert, Marvel G RENEWAL 3 RF H215/535754	28 Sep 2004			Hood, Maynard F	Yes	New Patient New PG		

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








Rapport

Internet
Fax
SMS

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Home Monitoring Service Center - Home Monitoring Demo System 1.6.0 09/03/05 - Microsoft Internet Explorer

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


PAC post IDM


DAI Belos-VRT

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Event report: Ventricular pacing impedance < 250 ohm or > 1500 ohm					
Standard report		Individual report		Complete report	
Save / print (PDF)		Help			
Event report: Ventricular pacing impedance < 250 ohm or > 1500 ohm					
Patient ID		Twiddler-Pace			
Implant / Implant SN		Belos VR-T / DEMO791545405			
Patient device SN		DEMO532545405			
Date of Cardio Report		02-Oct-2003 02:18			
		Since last follow-up day 11-Sep-2003		Since implantation	
DETECTION					
Episodes in VT1 zone			0	0	
Episodes in VT2 zone			0	0	
Episodes in VF zone			0	1	
THERAPY					
ATPs started			0	0	
ATPs successful			0	0	
Shocks started			0	1	
Shocks successful			0	1	
Shocks aborted			0	0	
30 J shock ineffective			0	0	
			Value	Measured on	
BATTERY					
Status			BOL	02-Oct-2003	
Voltage			6.18 V		
LEADS					
Ventricular pacing impedance			< 200 ohm	02-Oct-2003	
Shock impedance			31 ohm	26-Aug-2003	
IMPLANT STATUS SUMMARY					
Status		Ven. pac. imp. out of range			
Remark		Follow-up recommended			

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
Event report: Ventricular pacing impedance < 250 ohm or > 1500 ohm					Help
Standard report		Individual report		Complete report	Save / print (PDF)
Event report: Ventricular pacing impedance < 250 ohm or > 1500 ohm					
Patient ID			Twiddler-Pace		
Implant / Implant SN			Belos VR-T / DEM0791545405		
Patient device SN			DEMO532545405		
Date of Cardio Report			10-Oct-2003 02:17		
			Since last follow-up day 11-Sep-2003	Since implantation	
DETECTION					
Episodes in VT1 zone		0	0		
Episodes in VT2 zone		0	0		
Episodes in VF zone		0	1		
THERAPY					
ATPs started		0	0		
ATPs successful		0	0		
Shocks started		0	1		
Shocks successful		0	1		
Shocks aborted		0	0		
30 J shock ineffective		0	0		
		Value	Measured on		
BATTERY					
Status		BOL	10-Oct-2003		
Voltage		6.18 V			
LEADS					
Ventricular pacing impedance		> 3000 ohm	10-Oct-2003		
Shock impedance		31 ohm	26-Aug-2003		
IMPLANT STATUS SUMMARY					
Status			Ven. pac. imp. out of range		
Remark			Follow-up recommended		

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